

ECLIPSE

Overbed Table



ECLIPSE

Overbed Table

The Eclipse Overbed Table is a "trimmed-down" version of our popular Triumph overbed table. Engineered to fit a tight budget, while still maintaining a level of quality, the Eclipse overbed table features an attractive thermoformed top that is designed to contain large spills, as well as a sturdy powder-coated H base, as well as infinite height adjustment. Spring assisted lift mechanism ensures easy height adjustment from 28" to 45".

Features

- Thermoformed top with spill containment
- 16-Gauge tubular steel column and H-base
- Mechanical Spring mechanism

Options

- Vanity tray with mirror
- Wire storage basket
- Single wheel rubber casters
- Powder-coat or chrome finish



	Depth	Width	Height	Weight Capacity
Table	16"	30"	29"-45"	25 Lb.



Contact Us:
Email: Info@statmedicalcanada.com
Phone: 1-866-783-4076