

SPLIT TOP

Overbed Table



SPLIT TOP

Overbed Table

The Split Top Overbed Table provides extra functionality for both the patient and the caregiver. The secondary top can be slid out from either side of the upper top. The lower top is released by a tab and can be moved with one hand. When the tab is pressed, the two tops move in opposite directions.

Features

- Thermoformed top with spill containment
- Pneumatic cylinder
- Infinite height adjustment between 29.5" and 44.5"
- Oversized 2" dual-wheel casters
- Slide out top for extra surface area

Options

- Single wheel rubber casters
- Powder-coated or chrome plated 12-gauge steel column and base
- H- or U- Style base
- Column mounted on left or right side
- Vanity tray



	Depth	Width	Height	Weight Capacity
Table	18"	32"	29.5" – 44.5"	25 Lbs.



Contact Us:
Email: Info@statmedicalcanada.com
Phone: 1-866-783-4076