

TRANSCEND
Overbed Table



TRANSCEND

Overbed Table

The Transcend overbed table combines functionality with sleek lines and rugged durability. Its low-profile U-base is compatible with late-model beds, stretchers, mobile recliners and wheelchairs. Its pneumatic cylinder provides smooth, almost silent operation and permits infinite height adjustment within a height range of 30" to 45". The dimensions of the top are 18" x 32" x 3/4". A high-impact polystyrene shroud covers the base while concealing and protecting the twin-wheel rubber casters. Model ST180-PR boasts a highly durable poly-resin edge that provides a soft, comfortable surface. It's recessed for spill containment and molded to hold a cup and tissue box. Model ST180-TH features a seamless thermoformed top with spill containment.



Features

- Gas cylinder height adjustment mechanism
- 25 lb weight Capacity
- Infinite height positioning
- 29.5"-44.5" height range
- Anodized Aluminum column
- 14-gauge steel base covered with molded ABS shroud

Options

- Thermoform or poly-resin top
- Vanity tray with mirror
- Custom poly-resin edge colors
 - o Minimum 30 Unit order

	Depth	Width	Height	Weight Capacity
Table	18"	32"	30"-45"	25 lb.



Contact Us:
 Email: Info@statmedicalcanada.com
 Phone: 1-866-783-4076