

LIFEDOP

L250 Obstetrical and Vascular Doppler



STATMEDICAL

"we take care"

LIFEDOP

L250 Obstetrical and Vascular Doppler

Lifedop L250

All the things that you love in a Doppler come packed into this small package! The loud and clear sound is paralleled only by the large easy to read display for either the F.H.R. or the B.P.M... The easy to use power button and volume slider are conveniently located in intuitive locations, making it natural and familiar to healthcare professionals. You will also notice that this Doppler is suitable for obstetrical or vascular application with the large selection of probes available. Each probe is designed intelligently for the application, using broad beam technology for simplicity of use.



PROBE SELECTION

TRANS-VAGINAL



2MHZ Late term OBS



3MHZ Early term OBS



Water Proof Late term OBS



4MHZ Deep Vascular



5MHZ Mid Vascular



8MHZ Peripheral Vascular



Intra-operative surgical

Toll Free: 866.783.4076
www.statmedicalcanada.com