

PC-COMP

Geriatric Positioning Chair



PC-COMP

Geriatric Positioning Chair

The PC-Comp geriatric positioning chair requires few operating instructions and is fully adjustable to allow for maximum comfort and sore prevention. With its multiple configurations, this model alone appeals to a broad customer base.

Features

- Articulated headrest
- Seat cushion made of high resilience foam and visco gel memory foam
- Tilt-in-space
- Adjustable Seat width
- Adjustable seat height
- Adjustable backrest width

Options

- Catheter drainage bag hook
- 1 1/2" side bolsters
- Removable armrests
- PRE-CURE® PRO cushion
- Additional seat cushion
- Adjustable tension backrest
- 4" lateral trunk support
- 6" lateral trunk support
- Adjustable structural lateral supports
- Adjustable shoulder supports
- Gas cylinder assisted back angle adjustment 98 to 135 degrees
- Additional backrest cushion
- Opaque activity tray
- Activity tray with Lexan window
- Padded leatherette tray cover
- Hemiplegic armrest pad
- Padded leatherette footboard cover
- Lateral leg protectors
- Elevating footrest cushion
- 13" wheels with central brake
- IV pole
- 13"-wheel side brakes
- Lateral tray straps
- Wrap-around tray strap
- Universal perineal abductor kit with black fixing mat



	Depth	Depth Reclined	Width	Height	Weight Capacity
Chair	34.5"	65"	25.75"	53.25"	250 lbs
Seat	15.5-19.5"	N/A	16.25-19.25"	19-21"	N/A



Contact Us:
 Email: Info@statmedicalcanada.com
 Phone: 1-866-783-4076