

PC-HARMONY

Patient Positioning Chair



PC-HARMONY

Patient Positioning Chair

The Unique design of the PC-Harmony patient positioning chair gives it great versatility. In addition to clinical applications, it is also suitable as a geriatric recliner with its many possible adjustments. Its comfort is such that it is widely used in palliative care.

Features


- Removeable high-resilience foam cushion
- Adjustable seat height and depth
- Push back recline
- Push bar
- Two swivel wheels with brakes.

Options

- Left or right lateral tray
- Adjustable elbow supports
- Adjustable IV pole, rear installation
- Adjustable IV pole, front installation
- Adaptor for IV pole with 1/2" round tube
- Elevating footrest support
- Set of four swivel wheels with brakes
- Self-storing, removeable footrests
- Side bolster cushions
- Opaque activity tray
- Activity tray with Lexan window
- Lateral support bracket for storing activity tray
- PRE-CURE® PRO cushion allows the insertion of speciality pressure management gels or foams
- Inclined armrest for blood sample
- Contoured headrest
- Leatherette padded tray cover
- Catheter drainage bag hook
- Lateral tray straps
- Wraparound tray straps
- Universal perineal abductor kit with black fixing mat
- 4" lateral trunk supports with backrest canvas
- 6" lateral trunk supports with backrest canvas
- Cervical support



| | Depth | Depth Reclined | Width | Height | Weight Capacity |
|-------|----------|----------------|--------|------------|-----------------|
| Chair | 36.5" | 70" | 28.5" | 46.5" | 250 lbs |
| Seat | 17-18.5" | N/A | 20.25" | 17.5-20.5" | N/A |



Contact Us:

Email: Info@statmedicalcanada.com

Phone: 1-866-783-4076