

PC-OPT

Geriatric Positioning Chair



PC-OPT

Geriatric Positioning Chair



The PC-Opt Positioning Chair is an extremely versatile chair with a multitude of accessories to create a perfect chair for almost any requirement. Unique to this chair is the fold down head section that allows easy access to the upper head and neck area from the back of the chair. It also has drop arms that allow easy access for transfers or patient care.

Features

- Contoured polyurethane-covered seat for better load distribution
- Gas-cylinder assisted tilt
- Adjustable blocking ring to limit recline from 102 to 145°
- Backrest built on elastic bands



PC-OPT

Geriatric Positioning Chair

Options

- Lowered seat 2"
- Seat pan for specialized cushions
- Side bolsters 1"
- PRE-CURE® PRO cushion (Optional Seat pan rec)
- Adjustable lateral trunk support 4"
- Adjustable lateral trunk support 6"
- Adjustable shoulder support
- Cervical support
- Contoured headrest
- Adjustable tension backrest
- Adjustable structural lateral supports
- Adjustable elbow supports (pair)
- Additional contoured and articulated headrest
- Adjustable wings for upper section of backrest
- Additional upper backrest
- Opaque activity tray
- Activity tray with Lexan window
- Padded tray cover
- Hemiplegic armrest
- Padded footboard cover
- Lateral leg protector
- Footboard extension 2"
- Elevating footrest cushion
- Extra-large footboard
- Adjustable IV pole
- Catheter drainage bag hook
- Lateral tray straps
- Wrap-around tray strap
- Universal perineal abductor kit with black fixing mat



| | PC – Opt - 18 | PC – Opt - 20 | PC – Opt - 22 | PC – Opt - 24 |
|--------------------------|---------------|---------------|----------------|---------------|
| Depth | | | 38" | |
| Seat Depth | | | 15 1/2" to 19" | |
| Depth Reclined | | | 68" | |
| Width | 23 1/2" | 26 1/2" | | 29 1/4" |
| Width with mobility Pack | 27 1/2" | 32 1/2" | | 35 3/4" |
| Seat Width | 17 3/4" | 20 1/2" | 22" | 24" |
| Height | | | 50" | |
| Seat Height | | | 20 1/2" | |
| Weight Capacity | 300 lbs | | | 450 lbs |

Contact Us:
 Email: Info@statmedicalcanada.com
 Phone: 1-866-783-4076