

# CARSON

TREATMENT RECLINER



RECLINERS

# CARSON TREATMENT RECLINER

## Features

- Limited lifetime warranty
- Field replaceable upholstery covers
- Side lever allows both caregiver and patient to operate footrest
- Infinite position recline mechanism
- Ergonomic handles for transport

## Options

- Wood, solid surface, or polyurethane arm caps
- Fold down arm
- Central locking casters
- Trendelenburg
- Slide-out foot tray
- Flip-up side table
- Flip-swivel-lock side table
- Moisture barrier
- IV Pole
- Foley bag hook



	Depth	Width	Seat Width	Height	Seat Height	Weight Capacity
CRR123F	36"	29"	23"	48"	20"	500 lbs
CRR126F	36"	32"	26"	48"	20"	500 lbs
CRR123FT	36"	29"-32"	23"	48"	20"	500 lbs
CRR126FT	36"	32"-35"	26"	48"	20"	500 lbs

## Statement of Line



CRR123F



CRR126F



CR123FT



CRR126FT

