

COMPANION RECLINER



COMPANION RECLINER

Features

- Pneumatic Height adjustment
- 3-position recline
- Minimal footprint mobile base
- Flexible
- Locking casters

Options

- Wood, solid surface, or polyurethane arm caps
- Fixed height
- Contrasting fabrics



	Depth	Width	Seat Width	Height	Seat Height	Weight Capacity
RC108	34"	27"	23"	45"-48"	18"-21"	350 lbs
RC108-FH	34"	27"	23"	45"	18"	350 lbs

Statement of Line



RC108



RC108-FH

