

CONNECT ELECTRIC RECLINER



CONNECT ELECTRIC RECLINER

Features

- 12 year warranty
- Electric height adjustment
- Mobile transport base
- Lays completely flat
- Easy lateral transfers

Options

- Stirrups
- Folding side tablet
- Swivel side tablet
- Fold-down arm
- Rolling foot tray
- Soft grip push bar
- Ergotron arm
- Drainage bag holder
- Heated seat and massage
- Central locking casters
- Straight or flared back
- Trendelenburg
- Battery backup
- 3 button foot control
- IV pole
- Push bar
- Paper roll dispenser



| | Depth | Width | Seat Width | Height | Seat Height | Weight Capacity |
|------------|-------|-------|------------|---------|-------------|-----------------|
| CONNECTE29 | 36" | 30" | 24" | 46"-59" | 18"-31" | 650 lbs |

Statement of Line



CONNECTE29

