

INVERNESS

Patient Treatment Recliner



INVERNESS

Patient Treatment Recliner

Rarely does a recliner evoke the intense passion of patients and clinicians like the Inverness 24-hour Treatment Recliner. The Inverness is a superbly crafted treatment recliner with a patent pending suite of intuitive, practical and durable innovations to address patient and clinician comfort for a completely new clinical seating experience.

Features

- Low effort bilateral controls for easier patient use
- Supports clinician comfort with front locking brake tabs and back relieving side table heights
- Enhances facility efficiency with attractive, durable construction that is easy to clean
- Trendelenburg
- Infinite position backrest
- Swing arms
- 5" Casters with footplate
- Push handle

Options

- Heat and massage
- IV pole
- Accessory hook
- Pivot tablet
- Polyurethane arm caps
- Headrest pillow
- Lumbar pillow
- Side pillow
- Side tablet



	Depth	Depth Reclined	Width	Height	Weight Capacity
Chair	36"	75.5"	34.5"	51.75"	500 lb.
Seat	20.75"	N/A	26"	25.5"	N/A



Contact Us:
Email: Info@statmedicalcanada.com
Phone: 1-866-783-4076