

RC301

Patient Treatment Recliner



RC301

Patient Treatment Recliner

The RC301 Recliner is generously sized, deeper and wider for patients up to 300 lbs., with a steel frame reinforced to meet the needs of larger patients. Durable five-inch casters roll easily during transport and raise the chair to a comfortable height for staff. The RC301 has extra room and is a great solution for nocturnal therapies and is ideal for oncology, pre/post op, dialysis, infusion therapy, patient rooms, and other treatment needs.

Features


- Infinite position recline
- Pressure reduction and temperature smart foam in the seat and back
- Heavy gauge reinforced steel frame
- Removeable seat
- 5" Casters, with two brakes
- Roll out foot tray
- One fold away side tablet
- Push handle
- Wall hugger

Options

- Trendelenburg
- Break-away arm
- Additional fold-away side tablet
- Cupholder
- Heat and massage
- Black over-the-lap tray
- Swing away side tablet
- Head, Lumbar, and side pillow
- Blood draw arms support
- Push handle cover
- Head flap protector
- Arm covers
- Rear back flap
- Monitor mount
- Nocturnal sleep pad and guards
- Footrest extender
- IV pole
- Oxygen tank holder
- Bag hook
- Personal items storage bag
- Paper roll dispenser



	Depth	Depth Reclined	Width	Height	Weight Capacity
Chair	47"	74"	32"	48"	300 lbs
Seat	18"	N/A	22.5"	22.5"	N/A



Contact Us:
Email: Info@statmedicalcanada.com
Phone: 1-866-783-4076