



**STAT**MEDICAL  
“we take care”

# IMEXDOP CT+

## DURABILITY

A design that has lasted more than 20 years and still has the revered recognition as the best doppler. The heavy duty case and cable make it extremely tough.

## THEFT DETERRENT

The size of the Imexdop CT+ makes it less likely to be pocketed and stolen. It also comes standard with a docking station that the doppler calls home making it less likely to be left behind on a countertop by mistake.

## AMAZING SOUND

If you want to hear clearly, the Imexdop CT+ is still the loudest and clearest on the market today. The dual speaker system is an industry leader.

## LONG LIFE

We have customers that have been using the same doppler for over 20 years.

**Book a trial today**



Imexdop CT+





### **Technical Specifications**

Weight ..... 1.2 kg (2.7 lbs)  
Dimensions ..... 22.4x10.9x15.7 cm (8.8x4.3x6.3 in)  
Power Source ..... 110V or 220V recharger  
Battery..... Replaceable 12 Volt NiMH battery

### **Warranty**

*5 years parts and 1 year labor against manufacture defects. The waterproof probe has a 1 year parts and labor warranty against manufacture defects.*

### **Probes:**

2 MHz obstetric probe ..... T200  
3 MHz obstetric probe ..... T300  
5 MHz vascular probe ..... T500  
8 MHz vascular probe ..... T800  
2 MHz waterproof obstetric probe .....WP20

### **Options**

Roll Stand with Basket ..... ST3  
IV Pole/Roll Stand mount kit ..... V100