

RC653

Patient Treatment Recliner



RC653

Patient Treatment Recliner

The RC653 is our basic workhorse in the clinical recliner line. Non-porous blow molded fold down side trays are mounted on both sides. Each chair in the RC653 series features a standard Trendelenburg release with handle accessible from either side. It also has a built-in ergonomic headrest, full coverage uni-directional elastic band suspension system in the back frame that provides exceptional comfort and an elastic seat base. All these features combine to create the ultimate in comfort, durability, and value. For use in Dialysis, Oncology, Acute care, and other clinical environments.



Features

- Dual fold down side tablets
- Trendelenburg release on back with easy access handle
- Built in ergonomic headrest with headrest cover
- Removeable side panels for easy cleaning
- 4 Locking 5" Casters
- Total lock and directional lock casters
- 3 position recline

Options

- Bariatric Version available
- Steel or nylon casters
- Leg extensions
- Heat and massage
- IV pole
- Clear vinyl armrest protector
- Lumbar pillow
- Side cushion

	Depth	Depth Reclined	Width	Height	Weight Capacity
Chair	39" 44"	76.5"	32.25" 36.25"	49"	275 lbs 450 lbs
Seat	20"	N/A	21" 25"	21"	N/A



Contact Us:
Email: Info@statmedicalcanada.com
Phone: 1-866-783-4076