

SIERRA

Patient Treatment Recliner



SIERRA

Patient Treatment Recliner

The Sierra is the perfect solution for small areas with the tiny footprint and complete recliner capability. It brings incredible comfort and elegance into the room without crowding out the available space.

Features

- Minimal seams for ease of cleaning
- Two position push back recline
- Mechanisms are made of steel with a protective powder coat finish
- Arms, legs and seat are all connected with metal to metal connections for optimal structural integrity
- No exposed wood, and all upholstered pieces are double stitched for a more durable finish and pleasing aesthetic finish
- Premium locking casters



	Depth	Depth Reclined	Width	Height	Weight Capacity
Chair	28"	67.5"	26.25"	47.25"	300 lb.
Seat	19.5"	N/A	20"	18.25"	N/A

Options

- Push handle
- Moisture barrier fabric
- Fire retardant polyurethane foam
- Contrasting fabrics



Contact Us:
Email: Info@statmedicalcanada.com
Phone: 1-866-783-4076