

VERITY

Patient Recliner



VERITY

Patient Recliner

The Verity recliner is a perfect blend of quality and function with a very competitive price tag. The simplicity of this design allows you to maintain a premium experience at a low cost while still impressing the user of the chair.

The verity recliner offers healthcare-grade durability and exudes an inviting, homelike aesthetic that brings residential warmth to clinical environments.

Features

- Limited Lifetime warranty
- Field-replaceable upholstery covers
- 500lb. weight capacity
- 4 locking casters
- 3 position recline with push back mechanism
- Ergonomically positioned push handles for patient transport
- Stitching method reduces seams for increased infection resistance.

Options

- Moisture barrier
- Optional wood or polyurethane arm caps
- Contrasting fabrics.



	Depth	Depth Reclined	Width	Height	Weight Capacity
Chair	28"	66"	29 "	43.75"	500 lb.
Seat	22"	N/A	22"	18.5"	N/A



Contact Us:
Email: Info@statmedicalcanada.com
Phone: 1-866-783-4076