

# SC227

Sleep Chair

---



# SC227

## Sleep Chair



The SC227 Sleeper chair is the newest chair added to our lineup and it is getting a lot of attention. With the small profile it is a space saver while continuing to look beautiful and provide the sit/sleep functionality to any room.

### Features

- Space saving design
- 2 Positions, sit and sleep
- 3" durable locking casters
- High density foam
- Back push handle for mobility
- Front 2 casters locking

	Depth	Depth Reclined	Width	Height	Weight Capacity
Chair	36"	85"	36"	37"	350 lbs
Seat	21"	N/A	24"	16"	N/A

### Options

- Fire retardant foam
- Industrial casters
- 4 Locking casters
- Tablet arm



Contact Us:  
Email: [Info@statmedicalcanada.com](mailto:Info@statmedicalcanada.com)  
Phone: 1-866-783-4076