

T-GLIDE METAL

Rocker Glider Chair



T-GLIDE METAL


Rocker Glider Chair



The Thera-Glide Metal chair safely allows users to benefit from the therapeutic effects of a rocking chair while minimizing the risk of falls during transfer. With its patented auto-lock mechanism permitting safer transfers, the Thera-Glide chair is exceptionally comfortable and extremely secure due to its contoured cushions, extended armrests with palm grips, anti-tip stable base, and height adjustment.

Features

- Self-Locking rocking mechanism
- Stable base with anti-tip feet
- Height adjustable
- Lifetime frame warranty
- Rear wheels for mobility



T-GLIDE METAL

Rocker Glider Chair

Options

- Adjustable seat depth and backrest angle
15"-18", 95-115°
- Upholstered side panels
- Adjustable height kit for 21" to 24"
- Polyurethane armrests
- Gliding metal ottoman
- Fixed metal ottoman
- Opaque activity tray
- Activity tray with Lexan window
- Lateral tray straps
- Wrap-around tray straps
- Side storage bag
- 1" seat insert
- Universal perineal abductor kit with black fixing mat
- 4" Lateral trunk supports with backrest canvas
- 6" lateral trunk supports with backrest canvas
- Cervical support
- PRE-CURE® PRO seat cushion
- 3 Backrest styles
 - Contemporary
 - Classic
 - Comfort Plus



	T-GLIDE Metal Small	T-GLIDE Metal Regular	T-GLIDE Metal Large
Depth		32"	
Seat Depth		17 1/2"	
Width	22"	26"	29"
Seat Width	16 1/2"	20 1/2"	23 1/2"
Height		44 - 47"	
Seat Height		16" - 19"	
Weight Capacity	300 Lbs		450 Lbs

Contact Us:
 Email: Info@statmedicalcanada.com
 Phone: 1-866-783-4076