

T-GLIDE TRANSITIONAL Rocker Glider Chair



T-GLIDE TRANSITIONAL Rocker Glider Chair

The new Thera-Glide™ Transitional chair brings a modern design and sophisticated look to the Thera-Glide™ family. Its transitional style combines the traditional style's warmth with the contemporary style's modernism. The Thera-Glide™ Transitional chair safely allows users to benefit from the therapeutic effects of a rocking chair while minimizing the risks of falls during transfer. With its patented auto-lock mechanism permitting safer transfers, the Thera-Glide™ Chair is exceptionally comfortable and extremely secure due to its contoured cushions, extended armrests with palm grips, anti-tip stable base, and height adjustment.




Features

- Extended armrests with ergonomic palm grips
- Removeable upholstery for easy cleaning
- Automatic infinite position lock
- Glide motion lock-out
- Anti-tip side stabilizers
- Adjustable seat height
- Sinuous loop spring seat constructions
- Matching color undercarriage cover

Options

- Rear wheels
- Backrest switch delete
- Upholstered wood ottoman with footboard
- Contrasting Fabrics

| | Depth | Width | Height | Weight Capacity |
|-------|-------|---------|------------------|-----------------|
| Chair | 33" | 27 1/4" | 42 1/2 - 45 1/2" | 350 Lbs |
| Seat | 18" | 22" | 16-19" | N/A |



Contact Us:
Email: Info@statmedicalcanada.com
Phone: 1-866-783-4076