

T-GLIDE WOOD

Rocker Glider Chair



T-GLIDE WOOD

Rocker Glider Chair



The Thera-Glide Wood chair safely allows users to benefit from the therapeutic effects of a rocking chair while minimizing the risk of falls during transfer. With its patented auto-lock mechanism permitting safer transfers, the Thera-Glide chair is exceptionally comfortable and extremely secure due to its contoured cushions, extended armrests with palm grips, anti-tip stable base, and height adjustment.

Features

- Extended armrests with ergonomic palm grips
- Removable upholstery for easy cleaning
- Automatic infinite position lock
- Glide motion lock-out
- Anti-tip side stabilizers
- Adjustable seat height
- Matching color undercarriage cover



T-GLIDE WOOD

Rocker Glider Chair

Options

- Adjustable seat depth and backrest angle
15"-18", 95-115° (for Regular, standard on Small and large)
- Opaque activity tray
- Activity tray with Lexan window
- Lateral tray straps
- Wrap-around tray straps
- Rear wheels
- Side storage bag
- 1" Seat insert
- Seat wedge insert
- Universal perineal abductor kit with black fixing mat
- 4" lateral trunk supports with backrest canvas
- 6" lateral trunk supports with backrest canvas
- Undercarriage cover delete
- Backrest switch Delete
- Cervical support
- Upholstered wood ottoman with footboard
- PRE-CURE® PRO seat cushion
- Contrasting fabrics
- Contemporary, Classic, or Comfort Plus backrest



	T-GLIDE Wood Small	T-GLIDE Wood Regular	T-GLIDE Wood large
Depth		33"	
Seat Depth		17"	
Width	22 1/2"	26 1/2"	29 1/2"
Seat Width	16 1/2"	20 1/2"	23 1/2"
Height		44 - 47"	
Seat Height		16" - 19"	
Weight Capacity	275 Lbs		450 Lbs

Contact Us:
 Email: Info@statmedicalcanada.com
 Phone: 1-866-783-4076