

Phlebotomy Chairs



The Spirit

This basic model features two infinitely adjustable armrests, durable plastic resin seat, and non-marring leveling feet. A clean and simple design that stands up to years of use.



SPECIFICATIONS

	With Storage Drawer	No Storage Drawer
Overall Width	35.25" (89.5 cm)	29" (73.7 cm)
Seat Width	18" (45.7 cm)	18" (45.7 cm)
Arm Height Range	26.5" (67.3 cm) – 32.5" (82.6 cm)	26.5" (67.3 cm) – 32.5" (82.6 cm)
Seat Height	20.5" (52 cm)	20.5" (52 cm)
Seat Depth	15" (38 cm)	15" (38 cm)
Seat Back Height	16" (40.6 cm)	16" (40.6 cm)
Overall Height	36" (91.4 cm)	36" (91.4 cm)
Overall Length	20.5" (52 cm)	20.5" (52 cm)
Weight Capacity	300 lbs. (136 kg)	300 lbs. (136 kg)

The Solace

The Solace features two infinitely adjustable armrests, an upholstered seat and back, and non-marring leveling feet. The flip-up contoured urethane armrest provides comfort, easy access, and added security for the patient. Available with or without a storage cabinet.



SPECIFICATIONS

	With Storage Drawer	No Storage Drawer
Overall Width	35.25" (89.5 cm)	29" (73.7 cm)
Seat Width	18" (45.7 cm)	18" (45.7 cm)
Arm Height Range	27.75" (70.5 cm) – 33.75" (85.7 cm)	27.75" (70.5 cm) – 33.75" (85.7 cm)
Seat Height	22.5" (57.2 cm)	22.5" (57.2 cm)
Seat Depth	18.5" (47 cm)	18.5" (47 cm)
Seat Back Height	24.5" (62.2 cm)	24.5" (62.2 cm)
Overall Height	44.5" (113 cm)	44.5" (113 cm)
Overall Length	27" (68.6 cm)	27" (68.6 cm)
Weight Capacity	300 lbs. (136 kg)	300 lbs. (136 kg)

The Unity

The Unity features two infinitely adjustable armrests, an upholstered seat and back, and a 450 lb. weight capacity. The flip-up contoured urethane armrest provides comfort, easy access, and added safety for the patient. Available with or without a storage cabinet.



SPECIFICATIONS

	With Storage Drawer	No Storage Drawer
Overall Width	40.25" (102.2 cm)	34" (86.4 cm)
Seat Width	24.25" (61.6 cm)	24.25" (61.6 cm)
Arm Height Range	27.75" (70.5 cm) – 33.75" (85.7 cm)	27.75" (70.5 cm) – 33.75" (85.7 cm)
Seat Height	22.5" (57.2 cm)	22.5" (57.2 cm)
Seat Depth	18.5" (47 cm)	18.5" (47 cm)
Seat Back Height	16" (40.6 cm)	16" (40.6 cm)
Overall Height	36.5" (92.7 cm)	36.5" (92.7 cm)
Overall Length	25" (63.5 cm)	25" (63.5 cm)
Weight Capacity	450 lbs. (204 kg)	450 lbs. (204 kg)